

## Product datasheet for **MR222305A1V**

### Mouse Snrpd3 (NM\_026095) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Snrpd3 (NM\_026095) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Snrpd3  
**Synonyms:** 1700043E15Rik; 2310009E13Rik; AW046420; SMD3  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR222305 representing NM\_026095  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCTATTGGTGTGCCGATTAAGTCTTGACGAGGCTGAAGGCCACATAGTGACGTGTGAGACCAACA  
CCGGGGAAGTATACCGAGGGAAGCTCATTGAAGCAGAGGACAACATGAACTGTCAGATGTCCAACATCAC  
AGTCACATACAGAGATGGTCGAGTGGCAGCTGGAACAGGTATATACGTGGCAGGAAGATCCGATTT  
CTGATTTTGCCTGACATGCTGAAAAATGCACCCATGTTAAAGAGCATGAAAAATAAAACCAAGGCTCAG  
GGGCCGCAGAGGAAAAGCTGCTATCCTGAAGGCCAAGTGGCCGCACGAGGGAGAGGACGTGGGATGGG  
ACGTGGGAACATCTTCCAGAAGCGAAGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR222305 representing NM\_026095  
Red=Cloning site Green=Tags(s)

MSIGVPIKVLHEAEGHIVTCETNTGEVYRGKLEAEDNMNCQMSNITVYRDGRVAQLEQVYIRGRKIRF  
LILPDMLKNAPMLKSMKNKNQSGSAGRGKAAILKAQVAARGRGRGMGRGNIFQKRR

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_026095



[View online »](#)

<b>ORF Size:</b>	378 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_026095.4</a></u> , <u><a href="#">NP_080371.1</a></u>
<b>RefSeq Size:</b>	604 bp
<b>RefSeq ORF:</b>	381 bp
<b>Locus ID:</b>	67332
<b>UniProt ID:</b>	<u><a href="#">P62320</a></u>
<b>Cytogenetics:</b>	10 C1
<b>MW:</b>	14.4 kDa