

## Product datasheet for **MR222307L4V**

### Gimap5 (NM\_175035) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gimap5 (NM_175035) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gimap5
Synonyms:	D630024P16; E230026N22Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175035
ORF Size:	924 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222307).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_175035.5</a> , <a href="#">NP_778200.1</a>
RefSeq Size:	1943 bp
RefSeq ORF:	927 bp
Locus ID:	317757
UniProt ID:	<a href="#">Q8BWF2</a>
Cytogenetics:	6 B2.3
Gene Summary:	This gene encodes a protein belonging to the GTP-binding superfamily and to the immun-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. In humans, the IAN subfamily genes are located in a cluster at 7q36.1. [provided by RefSeq, Jul 2008]



[View online »](#)