

## Product datasheet for **MR222352L4V**

### Sorbs1 (NM\_001034962) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sorbs1 (NM_001034962) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sorbs1
Synonyms:	2310065E01Rik; 9530001P15Rik; CAP; mKIAA1296; Sh3d5; SH3P12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001034962
ORF Size:	2142 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222352).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001034962.1</a> , <a href="#">NP_001030134.1</a>
RefSeq Size:	5282 bp
RefSeq ORF:	2145 bp
Locus ID:	20411
UniProt ID:	<a href="#">Q62417</a>
Cytogenetics:	19 34.25 cM
Gene Summary:	Plays a role in tyrosine phosphorylation of CBL by linking CBL to the insulin receptor. Required for insulin-stimulated glucose transport. Involved in formation of actin stress fibers and focal adhesions.[UniProtKB/Swiss-Prot Function]



[View online »](#)