

Product datasheet for **MR222358A1V**

Mouse Tnfrsf12a (NM_001161746) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Tnfrsf12a (NM_001161746) AAV Particle
Tag: Myc-DDK
Symbol: Tnfrsf12a
Synonyms: AI255180; C87282; Fn14; HIP1; TWEAK-R; TweakR
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR222358 representing NM_001161746
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCTTCGGCTTGGCCGCGTCTCTGCCGAGATCCTCGTGTGGGATTCGGCTTGGTGTGATGCGCG
 CCGCGGCCGGGAGCAAGCACCAGGCGCCGAGCACCTCTGCCACTTCAGGCTACTGTGGCCATTCT
 GGGGGGCGCTCTTAGTCTGGTCTGTTTGGCGCTGGTTTCTAGTTTCTGCTGGAGAAGATGCCGC
 CGGAGAGAAAAGTTTACTACCCCATAGAGGAGACTGGTGGAGAGGGCTGCCAGGTGTGGCACTGATCC
 AG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR222358 representing NM_001161746
 Red=Cloning site Green=Tags(s)

MASAWPRSLPQILVLGFLVLMRAAAGEQAPGAAAPPAHFRLWPILGGALSLVLVLAIVSSFLVWRRCR
 RREKFTTPIETGGEGCPGVALIQ

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001161746
ORF Size: 282 bp


[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001161746.1</u> , <u>NP_001155218.1</u>
RefSeq Size:	875 bp
RefSeq ORF:	285 bp
Locus ID:	27279
UniProt ID:	<u>Q9CR75</u>
Cytogenetics:	17 A3.3
MW:	10.5 kDa