

Product datasheet for **MR222372L4V**

Kcne3 (NM_001190870) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kcne3 (NM_001190870) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kcne3
Synonyms:	2210017H05Rik; MiRP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001190870
ORF Size:	312 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222372).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001190870.1 , NP_001177799.1
RefSeq Size:	1333 bp
RefSeq ORF:	312 bp
Locus ID:	57442
UniProt ID:	Q9WTW2
Cytogenetics:	7 E2



[View online »](#)

Gene Summary:

Ancillary protein that assembles as a beta subunit with a voltage-gated potassium channel complex of pore-forming alpha subunits. Modulates the gating kinetics and enhances stability of the channel complex. Assembled with KCNB1 modulates the gating characteristics of the delayed rectifier voltage-dependent potassium channel KCNB1. Associated with KCNC4/Kv3.4 is proposed to form the subthreshold voltage-gated potassium channel in skeletal muscle and to establish the resting membrane potential (RMP) in muscle cells. Associated with KCNQ1/KCLQT1 may form the intestinal cAMP-stimulated potassium channel involved in chloride secretion that produces a current with nearly instantaneous activation with a linear current-voltage relationship.[UniProtKB/Swiss-Prot Function]