

## Product datasheet for **MR222373L4V**

### Dcun1d4 (NM\_178896) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dcun1d4 (NM_178896) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dcun1d4
Synonyms:	AI836376
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178896
ORF Size:	876 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222373).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178896.5</a> , <a href="#">NP_849227.1</a>
RefSeq Size:	4143 bp
RefSeq ORF:	879 bp
Locus ID:	100737
UniProt ID:	<a href="#">Q8CCA0</a>
Cytogenetics:	5 C3.3



[View online »](#)

**Gene Summary:**

Contributes to the neddylation of all cullins by transferring NEDD8 from N-terminally acetylated NEDD8-conjugating E2s enzyme to different cullin C-terminal domain-RBX complexes which are necessary for the activation of cullin-RING E3 ubiquitin ligases (CRLs).  
[UniProtKB/Swiss-Prot Function]