

Product datasheet for MR222377L4

Dnajc1 (NM_001190817) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dnajc1 (NM_001190817) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Dnajc1
Synonyms:	4733401K02Rik; AA960110; D230036H06Rik; Dnajl1; ERdj1; ERj1p; MTJ1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222377).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001190817
ORF Size:	1656 bp



[View online »](#)

