

Product datasheet for **MR222379L3V**

Kirrel3 (NM_026324) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kirrel3 (NM_026324) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kirrel3
Synonyms:	1500010O20Rik; 2900036G11Rik; mKIAA1867; NEPH2; SST4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026324
ORF Size:	2298 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222379).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026324.3 , NP_080600.1
RefSeq Size:	3701 bp
RefSeq ORF:	2301 bp
Locus ID:	67703
UniProt ID:	Q8BR86
Cytogenetics:	9 A4



[View online »](#)

Gene Summary:

Synaptic adhesion molecule required for the formation of target-specific synapses (PubMed:23637329, PubMed:26575286). Required for formation of target-specific synapses at hippocampal mossy fiber synapses. Required for formation of mossy fiber filopodia, the synaptic structures connecting dentate granule and GABA neurons. Probably acts as a homophilic adhesion molecule that promotes trans-cellular interactions and stabilize mossy fiber filopodia contact and subsequent synapse formation (PubMed:26575286). Required for the coalescence of vomeronasal sensory neuron axons (PubMed:23637329). May be involved in the hematopoietic supportive capacity of stroma cells; the secreted extracellular domain is directly responsible for supporting hematopoietic stem cells (PubMed:12665856).
[UniProtKB/Swiss-Prot Function]