

## Product datasheet for **MR222382L3V**

### **Kcne3 (NM\_001190871) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Kcne3 (NM_001190871) Mouse Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | Kcne3  |
| Synonyms:                 | 2210017H05Rik; MiRP2   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_001190871   |
| ORF Size:                 | 312 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR222382).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001190871.2</a> , <a href="#">NP_001177800.1</a>  |
| RefSeq Size:              | 1097 bp  |
| RefSeq ORF:               | 312 bp   |
| Locus ID:                 | 57442  |
| UniProt ID:               | <a href="#">Q9WTW2</a>   |
| Cytogenetics:             | 7 E2   |



[View online »](#)

**Gene Summary:**

Ancillary protein that assembles as a beta subunit with a voltage-gated potassium channel complex of pore-forming alpha subunits. Modulates the gating kinetics and enhances stability of the channel complex. Assembled with KCNB1 modulates the gating characteristics of the delayed rectifier voltage-dependent potassium channel KCNB1. Associated with KCNC4/Kv3.4 is proposed to form the subthreshold voltage-gated potassium channel in skeletal muscle and to establish the resting membrane potential (RMP) in muscle cells. Associated with KCNQ1/KCLQT1 may form the intestinal cAMP-stimulated potassium channel involved in chloride secretion that produces a current with nearly instantaneous activation with a linear current-voltage relationship.[UniProtKB/Swiss-Prot Function]