

Product datasheet for MR222384L4

G6b (Mpig6b) (NM_001033221) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	G6b (Mpig6b) (NM_001033221) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Mpig6b
Synonyms:	AU023871; C6orf25; G6b; Ng31
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222384).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001033221
ORF Size:	591 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001033221.3 , NP_001028393.1
RefSeq Size:	2252 bp
RefSeq ORF:	594 bp
Locus ID:	106722
Cytogenetics:	17 B1
Gene Summary:	Inhibitory receptor that acts as a critical regulator of hematopoietic lineage differentiation, megakaryocyte function and platelet production (PubMed:23112346). Inhibits platelet aggregation and activation by agonists such as ADP and collagen-related peptide (By similarity). This regulation of megakaryocyte function as well as platelet production and activation is done through the inhibition (via the 2 ITIM motifs) of the receptors CLEC1B and GP6:FcRgamma signaling (PubMed:23112346). Appears to operate in a calcium-independent manner (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR222384L4