

Product datasheet for **MR222467**

Tmem232 (NM_001008973) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem232 (NM_001008973) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tmem232
Synonyms:	E130009J12Rik; Gm940; Tes13-L; Tes13-S
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR222467 representing NM_001008973
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGGATCGCC

ATGAAAGCTTATAAGCCACCTGTGGTCAATAAGTTGGAGTAATATCTTCCACTTATCATGAAGAACTCC
TGAAATCTATTTTGAAGTTCAAACAGAAAGAGTCAAAAACCAAAGCCATCATTTTCAATTTCCAA
AGAATTCATATTGAGATTCAACACACAGATAATCCAGCAGAGGAGGAAGAACTGTTGGAACAAGCCAGG
AGACTCATCGTCAGGTCTAAGAGGAAGTTGGGCCTCAAGACACTGGGCTCTGGAAGCATGTGCATCTTC
CGACCGCATGGGCTGAGGTCACTATCTGGCTCAGTGCAAAGGGGAGATCCAAGTGAAGCTTTAAATAT
GCTCCACGCCCTTGGACCATGTCTCTTTGATCAGCATCAGCTTCTGCCCTGTTTTTCTCGCGGAA
TCCGTTTTATATCGACTCTGTTGCGATGCATTCATGAAGGGATTTTATATTTCAGTTGAGATAAACTGG
TGAAGATTGGCTATTTGATCTTCTACGGCTTTTGTATTCTTCTGCATGGTCACTTAGAAAGCTTTAA
ACAACATTTGCTTAGGCTGCAGCCATTTTATACGCACCTCATTCTCGAAACCGTCTACTACAAGTAT
CCGAACATCATCTCCAAGTGCAGTTCATCCTGAAAACATCGGAAATTATCTGTAAAAGAGAGCTTCATT
CTGAACCATTCGTGAAAAGCCCAGATGAACTGAGGATCCATATTCTGACCTGAATCATTTCGACGTGAA
TAAGAGGGGCTACGAAGTCAACCACCTGCTCTGGCACAGTGTGGCCGCTGGTCTCGTTCAGAACAAT
CGTCTCAGCTCACCGAAGTGTGGAACACCTCTCTTACAAAACGCAGCTTCAGACGAAATGCTGGT
TGGATTACGCCCTGGCTTTGATGGTCTGGGGAAAGCTGCTAAGCTGGACATGGCCTGCTGAAAACACT
GATGGACCTAGTGACAGATTTCTTGAGAACATCTGTCTGCTCAGAACCAGGAAGAGAAGTATAATATA
TATGATACTTCTGGGCCTCAGAGATTGATTCACATATAACAATTATTGCAGAAGTCTGTCTGTATG
CAGCTACCTCTGACCTGAGGAAAACAGCTTTGATTGGTTCTGTGCCTGTAAAAGTCCGCAGCAAGGTAT
TTCCCTCACAGACAAGTCTGAAGAGCTGCCAGAGTTGGATGGCGCAAGCATTTTAACTCTTCTGAAGTAT
TTCTCTCCAGAATCTCTGATAACTGTGAGAAGGTGATCTGGATTGGATACTATGGCATAGTGATAACC
TGGTGAAGATGTCATGGGAATTCAGGGTGAAGGATCAGGATGGACTTAGAAACATGATATGGCAAAC
ACTGAAAAGATAAAGGATTATGAGCAAGACCCACGGATCCGATGTGCATTAGTTATAGCCAGGCCGAA
TTAAATGGCCATCTGACCCTTCTGCACAAAGGCCACACCAAACCTCAGGAGAAGAAGTTTTCTCCAAGT
ACATCGGGTGGAGAATCGCTACTCTCTCCAGACTCTTTTCCCTCTCTCGATGTTGCCCTCCAAA
GACGCCAGTGAAGTGGACCTGCCAAGAAAACACACCATCAGAGAACGCCAGCCTGCAAAAAGAGAGTT
CTGCGCTTATCTTGAAGGATCATTCTTCTGTGGTAGAAGTTCCATGACTCCGATCCAATTTCTTCA
CCAAGGCAGATAAGAAATTGGAAGAAATCATTGATCATCATTGGCAGAAAGATATGGAGGCCGAAAGCG
AGAAGAGGAAGCATATAAGGCCAAAATCAGAAGGACAAGGAGGAAAAGGAGAAGATTCACTTTCAAGAA
ATTATGAAGCAAAGAGAACGGAAACTGAACAAACAAACCAACCCATGAAATAATTCTGTGAGAAAAG
AAAGTGGCTCAGAAAAGAAATGTGGGTTCTTCGAGCTGAAACCCAGCACCACCAATGCTAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR222467 representing NM_001008973
 Red=Cloning site Green=Tags(s)

```

MKAYKPPVVKFGVISSTYEELLKSFESSNRRKSQKPKPSFSISKEFILRFNHTDNPAAAAEELLEQAR
RLIVRSKRKLGKLTGSGKHVHLPTAWAEVIYLAQCKGEIQDEALNMLHASLDHVSFDHDQLPALFFLAE
SVLYRLCCDAFMKGYLYSVEIKLVKIGYIFLRLVFFFLHGHLESFKQHLLRLQPYLYALHFSEPSYYKY
PNIISNVQFILKTSEIICKRELHSEPFVESPEDETEDPYSDLNHLQLNKRGYEVNHLHWSVAWVSCVQNN
RPQLTEVLEHLLFYKTQLQTKCWLDSALALMVLGEAAKLDMACLKTMDLVDFLENILSAQNQEENYNI
YDTSWASEIVFTYTTIIAEVCLYAATSDLRKTALIGFCACKSPQQGISLTDKSEELPELDGASILTLLKY
FSSRISDNCEKVIWIGYIGIVYNLVKMSWELQGEQDGLRNMIWQTLQKIKDYEQDPRIRCALVIAQAE
LNGPSPDFCTKATPNSGEEVFSKYIGWRIATTLRSLFFPSLDVAPPKTPVEVDLPRKHTIRERQPAKRV
LRFILKDHSSVVEVSMTPYPNFFTKADKLEEIIDHHWQKDMARKREEEAYKAQNQKDKKEEKEIHFQE
IMKQREKRLNKQTKPYEIIILSEKESGSEKCGFFELKPSAPNAK
  
```

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9013_h08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001008973

ORF Size: 2025 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001008973.2](#), [NP_001008973.2](#)

RefSeq Size: 2608 bp

RefSeq ORF: 2028 bp


Locus ID: 381107

UniProt ID: [Q5K6N0](#)

Cytogenetics: 17 E1.1

MW: 78.7 kDa

Product images:



Circular map for MR222467