

Product datasheet for MR222467L4

Tmem232 (NM_001008973) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem232 (NM_001008973) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Tmem232
Synonyms:	E130009J12Rik; Gm940; Tes13-L; Tes13-S
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222467).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001008973
ORF Size:	2025 bp



[View online »](#)

