

## Product datasheet for **MR222476L4V**

### Spata24 (NM\_027733) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Spata24 (NM_027733) Mouse Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | Spata24  |
| Synonyms:                 | 2700012K08Rik; 4930583E11Rik; 5133400G04Rik; AU016220; TIPT; TIPT2   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_027733  |
| ORF Size:                 | 555 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR222476).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_027733.5</a> , <a href="#">NP_082009.3</a>  |
| RefSeq Size:              | 726 bp   |
| RefSeq ORF:               | 558 bp   |
| Locus ID:                 | 71242  |
| UniProt ID:               | <a href="#">Q6P926</a>   |
| Cytogenetics:             | 18 B2  |



[View online »](#)

**Gene Summary:**

Binds DNA with high affinity but does not bind to TATA boxes. Synergises with GMNN and TBP in activation of TATA box-containing promoters and with GMNN and TBPL1 in activation of the NF1 TATA-less promoter. May play a role in cytoplasm movement and removal during spermiogenesis.[UniProtKB/Swiss-Prot Function]