

Product datasheet for **MR222495L3V**

Sema3g (NM_001025379) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sema3g (NM_001025379) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sema3g
Synonyms:	AK129018
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001025379
ORF Size:	2340 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222495).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025379.1 , NP_001020550.1
RefSeq Size:	3686 bp
RefSeq ORF:	2343 bp
Locus ID:	218877
UniProt ID:	Q4LFA9
Cytogenetics:	14 B
Gene Summary:	Has chemorepulsive activities for sympathetic axons. Ligand of NRP2.[UniProtKB/Swiss-Prot Function]



[View online »](#)