

## Product datasheet for **MR222500L3V**

### Cep55 (NM\_028760) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cep55 (NM_028760) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cep55
Synonyms:	1200008O12Rik; 2700032M20Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028760
ORF Size:	1383 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222500).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028760.3</a>
RefSeq Size:	2735 bp
RefSeq ORF:	1386 bp
Locus ID:	74107
UniProt ID:	<a href="#">Q8BT07</a>
Cytogenetics:	19 C2



[View online »](#)

**Gene Summary:**

Plays a role in mitotic exit and cytokinesis. Recruits PDCD6IP and TSG101 to midbody during cytokinesis. Required for successful completion of cytokinesis. Not required for microtubule nucleation. Plays a role in the development of the brain and kidney.[UniProtKB/Swiss-Prot Function]