

Product datasheet for **MR222503L3V**

Cdk2ap1 (NM_013812) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdk2ap1 (NM_013812) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdk2ap1
Synonyms:	Apc10; Cdkap1; doc-1; Doc1; DORC1; p12; ST19
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013812
ORF Size:	345 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222503).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013812.2 , NP_038840.2
RefSeq Size:	1283 bp
RefSeq ORF:	345 bp
Locus ID:	13445
UniProt ID:	O35207
Cytogenetics:	5 F
Gene Summary:	specific inhibitor of the cell-cycle kinase CDK2.[UniProtKB/Swiss-Prot Function]



[View online »](#)