

Product datasheet for MR222506L3

Ubp1 (NM_001083319) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ubp1 (NM_001083319) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Ubp1
Synonyms:	Cp2b; LBP-1a; LBP-1b; NF2d9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222506).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001083319
ORF Size:	1620 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001083319.1
RefSeq Size:	3857 bp
RefSeq ORF:	1623 bp
Locus ID:	22221
UniProt ID:	Q811S7
Cytogenetics:	9 F3
Gene Summary:	Functions as a transcriptional activator in a promoter context-dependent manner. Involved in regulation of the alpha-globin gene in erythroid cells. Activation of the alpha-globin promoter in erythroid cells is via synergistic interaction with TFCP2. Functions as a trans-acting factor that regulates the domestic strain CYP2D9 gene through specific association with the regulatory element SDI-A1. Binding to SDI-A1 depends on the type of nucleotide at position 299; binding is abolished by a nucleotide substitution at this position. Modulates the placental expression of CYP11A1 (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR222506L3