

Product datasheet for **MR222507L3V**

Cst7 (NM_009977) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cst7 (NM_009977) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cst7
Synonyms:	Cmap
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009977
ORF Size:	501 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222507).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009977.3
RefSeq Size:	986 bp
RefSeq ORF:	435 bp
Locus ID:	13011
UniProt ID:	O89098
Cytogenetics:	2 G3
Gene Summary:	Inhibits papain and cathepsin L but with affinities lower than other cystatins. May play a role in immune regulation through inhibition of a unique target in the hematopoietic system. [UniProtKB/Swiss-Prot Function]



[View online »](#)