

### Product datasheet for MR222533L4

# Adam15 (NM\_009614) Mouse Tagged Lenti ORF Clone

#### **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** Adam15 (NM\_009614) Mouse Tagged Lenti ORF Clone

Tag: mGFP

Symbol: Adam15

Synonyms: AD56; MDC15; metar

Mammalian Cell Puromycin

Selection:

Vector:

pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this of

Sequence:

The ORF insert of this clone is exactly the same as(MR222533).

**Restriction Sites:** Sgfl-Mlul

Cloning Scheme:





<sup>\*</sup> The last codon before the Stop codon of the ORF

**ACCN:** NM\_009614

ORF Size: 2445 bp



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

#### Adam15 (NM\_009614) Mouse Tagged Lenti ORF Clone - MR222533L4

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 009614.2</u>

 RefSeq Size:
 2859 bp

 RefSeq ORF:
 2448 bp

 Locus ID:
 11490

 UniProt ID:
 088839


Cytogenetics: 3 39.07 cM

Gene Summary: This gene encodes a member of a disintegrin and metalloprotease (ADAM) family of endoproteases that play important roles in various biological processes including cell

signaling, adhesion and migration. This gene is prominently expressed in vascular cells, the endocardium, hypertrophic cells in developing bone, and specific areas of hippocampus and cerebellum. The encoded preproprotein undergoes proteolytic processing to generate a mature, functional protein. Mice lacking the encoded protein have increased bone mass resulting from osteoblast proliferation, and exhibit reduced neovascularization in a mouse model for retinopathy. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar processing. [provided by RefSeq, May 2016]



## **Product images:**



Circular map for MR222533L4