


Plasmid Map:



ACCN: NM_001177951

ORF Size: 2916 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001177951.1](#), [NP_001171422.1](#)

RefSeq Size: 3187 bp

RefSeq ORF: 2919 bp

Locus ID: 19893

Cytogenetics: X A1.1

Gene Summary: Could be a guanine-nucleotide releasing factor (By similarity). Plays a role in ciliogenesis (By similarity). Probably regulates cilia formation by regulating actin stress filaments and cell contractility (By similarity). May be involved in microtubule organization and regulation of transport in primary cilia (By similarity). Plays an important role in photoreceptor integrity. Isoform 5 may play a critical role in spermatogenesis and in intraflagellar transport processes. [UniProtKB/Swiss-Prot Function]