

Product datasheet for **MR222550L3V**

Tpbg (NM_001164792) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tpbg (NM_001164792) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tpbg
Synonyms:	5T4; AW495680; WAIF1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164792
ORF Size:	1278 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222550).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001164792.1 , NP_001158264.1
RefSeq Size:	3509 bp
RefSeq ORF:	1281 bp
Locus ID:	21983
UniProt ID:	Q9Z0L0
Cytogenetics:	9 E3.1
Gene Summary:	May function as an inhibitor of Wnt/beta-catenin signaling by indirectly interacting with LRP6 and blocking Wnt3a-dependent LRP6 internalization.[UniProtKB/Swiss-Prot Function]



[View online »](#)