

Product datasheet for **MR222593L3V**

Slitrk5 (NM_198865) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slitrk5 (NM_198865) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slitrk5
Synonyms:	2610019D03Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_198865
ORF Size:	2871 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222593).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198865.1 , NP_942565.1
RefSeq Size:	4829 bp
RefSeq ORF:	2874 bp
Locus ID:	75409
UniProt ID:	Q810B7
Cytogenetics:	14 E4
Gene Summary:	Suppresses neurite outgrowth.[UniProtKB/Swiss-Prot Function]



[View online »](#)