

Product datasheet for MR222594L2

Csf2 (NM_009969) Mouse Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Csf2 (NM_009969) Mouse Tagged Lenti ORF Clone

Tag: mGFP Symbol: Csf2

Synonyms: CSF; Csfgm; Gm-CSf; GMCSF; MGI-IGM

Mammalian Cell None

Selection:

Vector:pLenti-C-mGFP (PS100071)E. coli Selection:Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(MR222594).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_009969

ORF Size: 423 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Csf2 (NM_009969) Mouse Tagged Lenti ORF Clone - MR222594L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

P01587

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 009969.4</u>

RefSeq Size:1033 bpRefSeq ORF:426 bpLocus ID:12981

UniProt ID:

Cytogenetics: 11 32.13 cM


Gene Summary: Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from

various lineages, including granulocytes, macrophages, eosinophils and erythrocytes.


[UniProtKB/Swiss-Prot Function]



Product images:



Circular map for MR222594L2



Double digestion of MR222594L2 using Sgfl and Mlul