

Product datasheet for **MR222602L3V**

Dlg2 (NM_011807) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dlg2 (NM_011807) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dlg2
Synonyms:	A330103J02Rik; B230218P12Rik; B330007M19Rik; Dlg2; Gm1197; Gm21505; PSD93
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011807
ORF Size:	2556 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222602).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011807.3 , NP_035937.2
RefSeq Size:	7409 bp
RefSeq ORF:	2559 bp
Locus ID:	23859
UniProt ID:	Q91XM9
Cytogenetics:	7 51.07 cM



[View online »](#)

Gene Summary:

Required for perception of chronic pain through NMDA receptor signaling. Regulates surface expression of NMDA receptors in dorsal horn neurons of the spinal cord. Interacts with the cytoplasmic tail of NMDA receptor subunits as well as inward rectifying potassium channels. Involved in regulation of synaptic stability at cholinergic synapses. Part of the postsynaptic protein scaffold of excitatory synapses.[UniProtKB/Swiss-Prot Function]