

## Product datasheet for **MR222606L4V**

### Ctbp2 (NM\_009980) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctbp2 (NM_009980) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctbp2
Synonyms:	AA407280; D7Ertd45e; Gtrgeo6; Ribeye
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009980
ORF Size:	1335 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222606).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009980.4</a> , <a href="#">NP_034110.1</a>
RefSeq Size:	3047 bp
RefSeq ORF:	1338 bp
Locus ID:	13017
UniProt ID:	<a href="#">P56546</a>
Cytogenetics:	7 76.32 cM
Gene Summary:	Corepressor targeting diverse transcription regulators. Isoform 2 probably acts as a scaffold for specialized synapses (By similarity). Functions in brown adipose tissue (BAT) differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)