

## Product datasheet for **MR222639L4V**

### Fpr3 (NM\_008042) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Fpr3 (NM_008042) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fpr3
Synonyms:	ALX; Fpr-rs1; Fpr-s1; Fprl1; LXA4-R; Lxa4r
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008042
ORF Size:	1053 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222639).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008042.2</a> , <a href="#">NP_032068.2</a>
RefSeq Size:	1220 bp
RefSeq ORF:	1056 bp
Locus ID:	14294
UniProt ID:	<a href="#">O08790</a>
Cytogenetics:	17 A3.2
Gene Summary:	Low affinity receptor for N-formyl-methionyl peptides. Receptor for lipoxin A4. May have an olfactory function associated with the identification of pathogens or of pathogenic states. [UniProtKB/Swiss-Prot Function]



[View online »](#)