

## Product datasheet for **MR222642L3V**

### Tfe3 (NM\_001105196) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tfe3 (NM_001105196) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tfe3
Synonyms:	bHLHe33; F830016E06Rik; Tcfe3; Tfe-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001105196
ORF Size:	1611 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222642).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001105196.1</a> , <a href="#">NP_001098666.1</a>
RefSeq Size:	3188 bp
RefSeq ORF:	1614 bp
Locus ID:	209446
UniProt ID:	<a href="#">Q64092</a>
Cytogenetics:	X 3.5 cM



[View online »](#)

**Gene Summary:**

Transcription factor that specifically recognizes and binds E-box sequences (5'-CANNTG-3') (PubMed:16936731). Efficient DNA-binding requires dimerization with itself or with another MiT/TFE family member such as TFEB or MITF (PubMed:16936731). In association with TFEB, activates the expression of CD40L in T-cells, thereby playing a role in T-cell-dependent antibody responses in activated CD4(+) T-cells and thymus-dependent humoral immunity (PubMed:16936731). Specifically recognizes the MUE3 box, a subset of E-boxes, present in the immunoglobulin enhancer (By similarity). It also binds very well to a USF/MLTF site (By similarity). May regulate lysosomal positioning in response to nutrient deprivation by promoting the expression of PIP4P1 (PubMed:29146937).[UniProtKB/Swiss-Prot Function]