

Product datasheet for **MR222662L4V**

B3gnt5 (NM_001159407) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	B3gnt5 (NM_001159407) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	B3gnt5
Synonyms:	beta3GnT5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001159407
ORF Size:	1128 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222662).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001159407.1 , NP_001152879.1
RefSeq Size:	4872 bp
RefSeq ORF:	1131 bp
Locus ID:	108105
UniProt ID:	Q8BGY6
Cytogenetics:	16 A3



[View online »](#)

Gene Summary:

Beta-1,3-N-acetylglucosaminyltransferase that plays a key role in the synthesis of lacto- or neolacto-series carbohydrate chains on glycolipids, notably by participating in biosynthesis of HNK-1 and Lewis X carbohydrate structures. Has strong activity toward lactosylceramide (LacCer) and neolactotetraosylceramide (nLc(4)Cer; paragloboside), resulting in the synthesis of Lc(3)Cer and neolactopentaosylceramide (nLc(5)Cer), respectively. Plays a central role in regulating neolacto-series glycolipid synthesis during embryonic development.
[UniProtKB/Swiss-Prot Function]