

Product datasheet for **MR222699L4V**

Ptrh2 (NM_001098810) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ptrh2 (NM_001098810) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ptrh2
Synonyms:	A230072116Rik; Bit1; CGI-147
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001098810
ORF Size:	546 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222699).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001098810.2
RefSeq Size:	2942 bp
RefSeq ORF:	546 bp
Locus ID:	217057
UniProt ID:	Q8R2Y8
Cytogenetics:	11 C
Gene Summary:	The natural substrate for this enzyme may be peptidyl-tRNAs which drop off the ribosome during protein synthesis.[UniProtKB/Swiss-Prot Function]



[View online »](#)