

Product datasheet for **MR222703L3V**

Sema3b (NM_009153) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sema3b (NM_009153) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sema3b
Synonyms:	LUCA; Se; sem; SemA; sema5; Sema; semaV
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009153
ORF Size:	2247 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222703).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009153.2 , NP_033179.2
RefSeq Size:	3115 bp
RefSeq ORF:	2250 bp
Locus ID:	20347
Cytogenetics:	9 58.31 cM



[View online »](#)

Gene Summary:

This gene encodes a secreted protein that belongs to the class 3 semaphorin/collapsin family. Members of this family play a role in growth cone guidance during neurogenesis. The encoded protein inhibits axonal extension. This protein is thought to be an osteoblast protein that regulates bone mass and affects skeletal homeostasis. A similar gene in humans functions as a tumor suppressor gene. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]