

Product datasheet for **MR222715L3V**

Ltb (NM_008518) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ltb (NM_008518) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ltb
Synonyms:	AI662801; LTbeta; p33; Tnfc; Tnfsf3; Tnlg1c
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008518
ORF Size:	918 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222715).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008518.2 , NP_032544.1
RefSeq Size:	1122 bp
RefSeq ORF:	921 bp
Locus ID:	16994
UniProt ID:	P41155
Cytogenetics:	17 18.59 cM
Gene Summary:	Cytokine that binds to LTBR/TNFRSF3. May play a specific role in immune response regulation. Provides the membrane anchor for the attachment of the heterotrimeric complex to the cell surface.[UniProtKB/Swiss-Prot Function]



[View online »](#)