

Product datasheet for MR222717

Lamc2 (NM_008485) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lamc2 (NM_008485) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lamc2
Synonyms:	AA589349
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR222717 representing NM_008485 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGCGCTCTGGCTCAGCTGCTGCCTTGGTGTGCGCTCCTGCTGCCCGCAGCCAGGCCACCTCCA
GGAGGGAAGTCTGTGATTGCAATGGGAAGTCCAGGCAATGTGTCTTTGATCAGGAGCTCCATCGACAGAC
AGGCAGCGGGTCCGTTGCCTCACTGCAATGACAATACAGCGGGGTTCACTGCGAGCGGTGCAGGGAG
GGGTTTTACCGACATCGAGACAGAGACCGCTGCCTACCCTGCAACTGCCACTCAAAGGGTCCCTCAGTG
CTGGATGTGACAACTCTGGACAATGCAGGTGTAAGCCAGGTGTGACAGGACAAAGATGTGACCGATGTCA
GCCAGGCTTCATATGCTCACCGATGCTGGATGCACCCGAGACCAGGGGCAACTAGATTCCAAGTGTGAC
TGTGACCCAGCTGGCATCTCTGGACCCTGTGATTCTGGCCGATGTGTCTGCAAACCAGCCGTCACTGGAG
AGCGCTGTGATAGGTGCCGACCAGGCTACTATCATCTGGACCGGGCAAACCTGAGGGCTGTACCCAGTG
TTTCTGCTATGGGCATTAGCCAGCTGCCACGCCCTCTGCCGACTTCAGTGTCCACAAAATCACTTCAACT
TTCAGTCAGGATGTGGATGGTTGGAAGGCGGTTACAGAAAACGGGGCACCTGCAAAACTCCACTGGTCAC
AGCGCCATCGGGACGTGTTAGTTCTGCCGAAGATCAGACCCGCTATTTCTGGCCCTGCCAAATT
CCTCGGTAACCAAGTGAAGTACGGGACAGAGCTGTCTTTGACTACCGGTGGACAGAGGAGGTAGA
CAGCCGCTGCCTACGATGTGATCCTGGAAGGTGCTGGTCTACAGATCAGAGCTCCTGTGATGGCTCCAG
GCAAGACACTTCCCTGTGGGATCACAAGACTTACACATTACAGACTGAATGAACATCCAAGCAGTCACTG
GAGTCCCCAGCTGAGTTATTTGCAATATCGAAGGTTACTGCGGAACCTCACAGCCCTCCTGATCCGAGCT
ACGTACGGAGAATATAGTACAGGGTACATTGATAACGTGACCCTGGTTTCAGCCCGCCCTGTCTCTGGAG
CCCCAGCCCTTGGGTTGAACGTTGTGTATGCCCTGCTGGGTACAAGGGACAATTCTGCCAGGAATGTGC
TTCTGGTTACAAAAGAGATTCGGCAAGATTGGGGCCTTTTGGCGCTGTGTTCCCTGTAACCTGCCAAGGG
GGAGGGGCTGTGATCCAGACACGGGAGATTGCTACTCGGGGACGAGAATCCTGACATTGAGTGTGCTG
ACTGTCCCATCGTTTTACAATGACCCACATGACCCCGCAGCTGCAAGCCATGTCCCTGTCAATGG
GTTGAGCTGTTGAGTGCCTGAGACAGAGGAGTGGTGTGTAACAACCTGTCCCTGGGGTACAGGT



[View online »](#)

GCCCCGTGTGAGCTCTGTGCTGATGGCTTCTTTGGGGATCCCTTTGGGGAACGTGGCCCAGTGAGGCCTT
GTCAACGCTGCCAATGCAACAACAACGTGGACCCCAATGCCTCTGGGAACTGTGACCAGTTGACAGGCAG
ATGCTTGAATGTATCTACAACACGGCCGGTGTCTACTGTGACCAGTGCAAAGCAGGTTACTTTGGAGAC
CCATTGGCTCCCAACCCAGCAGACAAGTGTGAGCTTGAACCTGCAGCCCATGGGCTCGGAGCCTGGAG
AGTGTGAGGTGATGGCAGCTGTGTTTGAAGCCAGGCTTTGGCGGCCTCAACTGTGATCACGCAGCCCT
AACCAGTTGCTCTGCTTGTACAATCAAGTGAAGATTGAGTGGACCAGTTTACCAGCAGCTCCAGAGC
CTGGAGGCCCTGGTTTTCAAAGGCTCAGGGTGGTGGTGGTACAGTCCCCAGTGGGAGCTGGAAG
GCAGGATGCAGCAGGCTGAGCAGGCCCTTCGAGACATTCTAAGAGAAGCTCAGATTTAGAAAGGGCAAT
GAGAGCCCTTGGCTTCCGGCTGGCCAAGGCAAGGAGCCAAAGAGAACGACTACAAGACCCGCTGGATGAC
CTCAAGATGACTGCAGAAAGGATCCGGGCCCTGGGCAGTCAGCATCAGAACAGAGTTGAGGATACGAGCA
GACTCATCTCTCAGATGCGCCTGAGTCTGGCAGGAAGCGAAGCTTCCCTGCAAAACACTAATATCCATTC
TTCTGAGCACTACGTGGGCGCAATGGATTTAAAAGTCTGGCTCAGGAGGCTACAAGATTGGCAGACAGC
CACGCTGAGTCAGCTAACGCAATGAAGCAACTAGCAAGGAAAAGTGGAGTACTCCAAACAGCACTTT
CATTGGCCCGCAAGGCTCTGAGTGGAGGAGGCGGAAGTGGCGTCTTGGACAGCTCCGTGGTACAAGTCT
TATGGGAAAATTAGAGAAAACCAAGTCCCTGAGCCAGCAGCTGTCAAGGGAGGGCACCAAGCCGACATT
GAAGCTGATAGGTCGTATCAGCACAGTCTCCGCCCTCTGGATTCTGCCTCTCAGTTCAGGGAGTCAAGT
ATCTGTCTTTCAGGTGGAAGCAAGAGGATCAGACAAAAGGCTGATTCTCTCTCAAACCTGGTGACCAG
ACAAACGGATGCATTCACGCGTGTGCGAAAACAATCTGGGAACTGGGAAAAAGAAACACGGCAGCTTTTA
CAGACTGGAAAGGATAGGAGACAGACTTCAGATCAGCTGCTTCCCGTGCCAACCTTGCTAAAAACAGAG
CCCAAGAAGCGCTAAGTATGGGCAATGCCACTTTTTATGAAGTTGAGAACATCCTGAAGAACCTCCGAGA
GTTTGATCTGCAGTTGAAGACAGAAAAGCAGAGGCTGAAGAGGCCATGAAGAGACTCTCCTCTATTAGC
CAGAAGTTGCAGATGCCAGTGACAAGACCCAGCAAGCAGAAACGGCCCTGGGGAGCGCCACTGCCGACA
CCCAACGGGCAAAGAACGCAGCTAGGGAGGCCCTGGAGATCAGCAGCGAGATAGAGCAGGAGATAGGGAG
TCTGAACTTGAAGCTAATGTGACAGCAGATGGGGCCTTGGCCATGGAGAAGGGACTGGCCACTCTGAAG
AGCGAGATGAGAGAGATGGAGGGGAGCTGGCCAGAAAGGAGCTGGAGTTTGACACGGATAGGACACGG
TGCAGCTGGTGATTACTGAAGCCAGCAAGCTGATGCCAGAGCCACGAGTGCCGGAGTTACCATCCAAGA
CACACTCAACACATTGGACGGCATCCTACACCTCATAGACCAGCCTGGCAGTGTGGATGAAGAAGGGATG
ATGCTATTAGAACAAGGGCTTTTCCAAGCCAAGACCCAGATCAACAGTCGACTTCGGCCCTTGATGTCTG
ACCTGGAGGAGAGGGTGGTGGCAGAGGAACCACTCCATCTGCTGGAGACCAGCATAGATGGAATTCT
AGCTGATGTGAAGAACCTGGAGAACATTCGAGACAACCTGCCCCAGGCTGCTACAATACCAAGCTCTT
GAGCAACAG

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR222717 representing NM_008485
 Red=Cloning site Green=Tags(s)


MPALWLSCLGLVALLLPAAQATSRREVDCNGKSRQCVDQELHRQTGSGFRCLNCNDNTAGVHCERCRE
 GFYRHRDRDRCLPCNCHSKGSLSAGCDNSGQCRCCKPGVTGQRCDCRCQPGFHM LTDAGCTRDQGLDSKCD
 CDPAGISGPCDSGRVCVCKPAVTGERCDRCRPGYYHLDRANPEGCTQCFCYGHASCHASADF SVHKITST
 F SQD V D G W K A V Q R N G A P A K L H W S Q R H R D V F S S A R R S D P V Y F V A P A K F L G N Q Q V S Y G Q S L S F D Y R V D R G G R
 Q P S A Y D V I L E G A G L Q I R A P L M A P G K T L P C G I T K T Y T F R L N E H P S S H W S P Q L S Y F E Y R R L L R N L T A L L I R A
 T Y G E Y S T G Y I D N V T L V S A R P V S G A P A P W V E R C V P A G Y K Q F C Q E C A S G Y K R D S A R L G P F G A C V P C N C Q G
 G G A C D P D T G D C Y S G D E N P D I E C A D C P I G F Y N D P H D P R S C K P C P C H N G F S C S V M P E T E E V V C N N C P P G V T G
 A R C E L C A D G F F G D P F G E R G P V R P C Q R C Q C N N N V D P N A S G N C D Q L T G R C L K C I Y N T A G V Y C D C K A G Y F G D
 P L A P N P A D K C R A C N C S P M G S E P G E C R G D G S C V C K P G F G G L N C D H A A L T S C P A C Y N Q V K I Q M D Q F T Q Q L Q S
 L E A L V S K A Q G G G G G T V P S G E L E G R M Q Q A E Q A L R D I L R E A Q I S E G A M R A L G F R L A K A R S Q E N D Y K T R L D D
 L K M T A E R I R A L G S Q H Q N R V Q D T S R L I S Q M R L S L A G S E A S L Q N T N I H S S E H Y V G P N G F K S L A Q E A T R L A D S
 H A E S A N A M Q L A R E T E D Y S K Q A L S L A R K A L S G G G G S G V L D S S V V Q G L M G K L E K T K S L S Q Q L S R E G T Q A D I
 E A D R S Y Q H S L R L L D S A S Q L Q G V S D L S F Q V E A K R I R Q K A D S L S N L V T R Q T D A F T R V R N N L G N W E K E T R Q L L
 Q T G K D R R Q T S D Q L L S R A N L A K N R A Q E A L S M G N A T F Y E V E N I L K N L R E F D L Q V E D R K A E A E E A M K R L S S I S
 Q K V A D A S D K T Q Q A E T A L G S A T A D T Q R A K N A A R E A L E I S S E I E Q E I G S L N L E A N V T A D G A L A M E K G L A T L K
 S E M R E M E G E L A R K E L E F D T D K D T V Q L V I T E A Q Q A D A R A T S A G V T I Q D T L N T L D G I L H L I D Q P G S V D E E G M
 M L L E Q G L F Q A K T Q I N S R L R P L M S D L E E R V R R Q R N H L H L L E T S I D G I L A D V K N L E N I R D N L P P G C Y N T Q A L
 E Q Q

SGPTRRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9007_d11.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_008485

ORF Size: 3579 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008485.3](#), [NP_032511.3](#)

RefSeq Size: 5182 bp


RefSeq ORF: 3582 bp

Locus ID: 16782

Cytogenetics: 1 65.26 cM

MW: 130.8 kDa

Product images:



Circular map for MR222717