

Product datasheet for **MR222768L4V**

Tmem266 (NM_172923) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem266 (NM_172923) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem266
Synonyms:	9630029F15; A1118078; HVRP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172923
ORF Size:	1614 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222768).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172923.3
RefSeq Size:	2407 bp
RefSeq ORF:	1617 bp
Locus ID:	244886
UniProt ID:	Q8BZB3
Cytogenetics:	9 B



[View online »](#)

Gene Summary:

Voltage-sensor protein present on the post-synaptic side of glutamatergic mossy fibers and granule cells in the cerebellum. Despite the presence of a voltage-sensor segment, does not form a functional ion channel and its precise role remains unclear. Undergoes both rapid and slow structural rearrangements in response to changes in voltage. Contains a zinc-binding site that can regulate the slow conformational transition.[UniProtKB/Swiss-Prot Function]