

Product datasheet for **MR222790L4V**

Prss53 (NM_001081268) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Prss53 (NM_001081268) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Prss53
Synonyms:	mFLJ00289; Pol3s
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001081268
ORF Size:	1656 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222790).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001081268.1 , NP_001074737.1
RefSeq Size:	2531 bp
RefSeq ORF:	1659 bp
Locus ID:	330657
UniProt ID:	Q571E5
Cytogenetics:	7 F3
Gene Summary:	In vitro can degrade the fibrinogen alpha chain of as well as pro-urokinase-type plasminogen activator.[UniProtKB/Swiss-Prot Function]



[View online »](#)