

Product datasheet for MR222799L3

S1pr4 (NM_010102) Mouse Tagged Lenti ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | S1pr4 (NM_010102) Mouse Tagged Lenti ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | S1pr4 |
| Synonyms: | Edg6; lp; Lpc1; S1; S1p4 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| E. coli Selection: | Chloramphenicol (34 ug/mL) |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR222799). |
| Restriction Sites: | SgfI-MluI |
| Cloning Scheme: | |

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


| | |
|-----------|-----------|
| ACCN: | NM_010102 |
| ORF Size: | 1158 bp |



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_010102.2 , NP_034232.1 |
| RefSeq Size: | 2393 bp |
| RefSeq ORF: | 1161 bp |
| Locus ID: | 13611 |
| UniProt ID: | Q9Z0L1 |
| Cytogenetics: | 10 C1 |
| Gene Summary: | This gene encodes a member of the G-protein coupled receptor 1 family. The encoded protein is a receptor for the lysophospholipid sphingosine 1-phosphate. The encoded protein is one of the predominant sphingosine 1-phosphate G-protein coupled receptors expressed in T cells, and it plays a role in suppression of T cell proliferation, generation of cytokines, and potentially other aspects of immune activation. [provided by RefSeq, Feb 2010] |

Product images:



Circular map for MR222799L3