

## Product datasheet for **MR222809L3V**

### Slc2a9 (NM\_001102414) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc2a9 (NM_001102414) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc2a9
Synonyms:	Glut9; SLC2a9A; SLC2A9B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001102414
ORF Size:	1617 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222809).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001102414.1</a> , <a href="#">NP_001095884.1</a>
RefSeq Size:	3602 bp
RefSeq ORF:	1617 bp
Locus ID:	117591
UniProt ID:	<a href="#">Q3T9X0</a>
Cytogenetics:	5 B3
Gene Summary:	Urate transporter, which may play a role in the urate reabsorption by proximal tubules (PubMed:19587147, PubMed:25100214). Does not transport glucose, fructose or galactose (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)