

Product datasheet for **MR222810L4V**

Sox17 (NM_011441) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sox17 (NM_011441) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sox17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011441
ORF Size:	1257 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222810).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011441.4
RefSeq Size:	3130 bp
RefSeq ORF:	1260 bp
Locus ID:	20671
UniProt ID:	Q61473
Cytogenetics:	1 1.65 cM



[View online »](#)

Gene Summary:

This gene encodes a member of the Sox (Sry-related high mobility group box) family of transcription factors involved in the regulation of embryonic development. The encoded protein plays a role in the determination of cell fate and in maintaining cell identity. This gene regulates tumor angiogenesis and tumor progression. Mutations in the human gene are associated with vesicoureteral reflux, characterized by the backward flow of urine from the bladder into the ureters or the kidney. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]