

## Product datasheet for **MR222814L3V**

### Plaur (NM\_011113) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Plaur (NM_011113) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Plaur
Synonyms:	Cd87; u-PAR; uPAR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011113
ORF Size:	981 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222814).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011113.3</a> , <a href="#">NP_035243.1</a>
RefSeq Size:	1396 bp
RefSeq ORF:	984 bp
Locus ID:	18793
UniProt ID:	<a href="#">P35456</a>
Cytogenetics:	7 A3
Gene Summary:	Acts as a receptor for urokinase plasminogen activator. Plays a role in localizing and promoting plasmin formation. Mediates the proteolysis-independent signal transduction activation effects of U-PA.[UniProtKB/Swiss-Prot Function]



[View online »](#)