

Product datasheet for **MR222834L4V**

Fiz1 (NM_011813) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fiz1 (NM_011813) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fiz1
Synonyms:	AI790204; mFLJ00416
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011813
ORF Size:	1500 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222834).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011813.3 , NP_035943.3
RefSeq Size:	2656 bp
RefSeq ORF:	1503 bp
Locus ID:	23877
UniProt ID:	Q9WTJ4
Cytogenetics:	7 A1
Gene Summary:	May be a transcriptional repressor of NRL function in photoreceptors. Does not repress CRX-mediated transactivation (By similarity).[UniProtKB/Swiss-Prot Function]


[View online »](#)