

Product datasheet for **MR222836L3V**

Krt18 (NM_010664) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Krt18 (NM_010664) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Krt18
Synonyms:	CK18; K18; Krt1-18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010664
ORF Size:	1269 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222836).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010664.2 , NP_034794.2
RefSeq Size:	1400 bp
RefSeq ORF:	1272 bp
Locus ID:	16668
UniProt ID:	P05784
Cytogenetics:	15 57.22 cM


[View online »](#)

Gene Summary:

When phosphorylated, plays a role in filament reorganization. Involved in the delivery of mutated CFTR to the plasma membrane. Involved in the uptake of thrombin-antithrombin complexes by hepatic cells (By similarity). Together with KRT8, is involved in interleukin-6 (IL-6)-mediated barrier protection.[UniProtKB/Swiss-Prot Function]