

## Product datasheet for **MR222839L4V**

### Xpa (NM\_011728) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Xpa (NM_011728) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Xpa
Synonyms:	AI573865; Xpac
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011728
ORF Size:	816 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222839).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011728.2</a> , <a href="#">NP_035858.2</a>
RefSeq Size:	972 bp
RefSeq ORF:	819 bp
Locus ID:	22590
UniProt ID:	<a href="#">Q64267</a>
Cytogenetics:	4 24.49 cM



[View online »](#)

**Gene Summary:**

Involved in DNA excision repair. Initiates repair by binding to damaged sites with various affinities, depending on the photoproduct and the transcriptional state of the region. Required for UV-induced CHEK1 phosphorylation and the recruitment of CEP164 to cyclobutane pyrimidine dimers (CPD), sites of DNA damage after UV irradiation (By similarity).[UniProtKB/Swiss-Prot Function]