

Product datasheet for **MR222851**

Abcd3 (NM_008991) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Abcd3 (NM_008991) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Abcd3
Synonyms:	AI313901; AU018866; AW146054; PMP; PMP68; PMP70; Pxm; Pxmp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR222851 representing NM_008991
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGCCCTTCAGCAAGTACTTGACGGCACGGAAACACCTCGCTGGCGGGGGCCGCTTCTGCTGCTCT
GCCTACTCCATAAGCGGGTTCGCGCCCTCGGCCTGCACGGTAAGAAAAGTGGAAAACCGCCATTACAGAA
TAATGAGAAAAGGAAAAGAGAGCGAGCTGTGGTGGACAAAGTGTTCATCAAGGCTCTCACAGATC
CTAAAAATTATGGTTCCTAGAACATTTTGTAAAGAGACAGGCTACTTGTACTTATTGCGGTTATGCTGG
TATCTCGAACATACTGTGATGTTGGATGATTGAGAATGGGACGCTCATTGAGAGTGGCATCATTGGTGC
CAGCAGCAAAGATTTCAAGAGATACTTATTCAACTTCATCGTGCATGCCTCTTATCTCTCTGGTTAAT
AACTTCTGAAGTATGGGTTAAATGAGCTCAAAGTGTGCTTCGAGTCAGGCTCACACGGTACCTCTACG
AGGAGTATCTCCAAGCTTTCACCTATTAAAAATGGGCAACCTGGATAACAGAATAGCAAATCCAGACCA
GCTGCTTACACAGGATGTAGAAAAATTTGTAACAGTGTAGTCGATCTGTATTCAAATCTTAGTAAGCCA
TTTTTAGACATAGTTTTATATATTTTCAAGTTAACAAGTGAATTTGGAGCTCAGGCCCCGCAAGCATGA
TGGCCTACTTGCTGTTTCTGGGCTATTCTAACTCGACTCAGAAGACCCATTGGTAAGATGACAATTAT
GGAGCAGAAGTACGAAGGAGAATACAGATACGTGAATTCACGGCTTATCACTAATAGTGAAGAAATTGCC
TTTTACAATGGGAATAACCGGGAGAAGCAGACAATCCACTCAGTCTTCGAAAACCTGGTGGAACACCTAC
ATAATTTCAATTTCTTTCGGTTTTCTATGGGTTTCATTGATAGCATCATTGCCAAATATGTTGCCACTGT
CGTCGGGTACCTGGTGTGAGTCGCCCCCTCTAGATCTGGCACACCCTCGCCACCTTACAGCACACAC
TCAGAGCTGCTGGAGGATTACTACAAAGTGAAGGATGCTTTTGAGAATGTCTCAAGCTTTGGGTCGTA
TAGTTTTGGCTGGGCGTGAATGACTAGATTGGCTGGTTTTACGGCTCGGATTACAGAATTAATGCAAGT
ACTAAAGGACTTAAATCATGGCAGATATGAACGTACAATGGTGTACAACAGGAGAAGGGTATTGAAGGA
GCACAAGCTAGTCCCTTGGTCCCTGGTGTGGAGAAATCATCAATACAGACAACATTATAAAGTTTGATC
ATGTTCTTTAGCAACACCAAAATGGTGACATCTTGATCCAAGACCTTAGTTTTGAAGTTCGATCTGGGGC
CAACGTTCTGATTTGGTCCAAACGGCTGTGGAAAGAGCTCCCTCTCCGTGTTCTTGGTGAGTTATGG
CCTTTTTTGGAGGGCGGCTTACTAAACCTGAGAGAGGAAAGTTATTTTATGTTCTCAGCGACCCTATA
TGACCCTGGGAACACTGAGAGACCAAGTCATATCCAGATGGAAAGGAAGATCAGAAGAAGAGGGGGAT
CTCTGACCAGGTGCTGAAGGAGTACTTGGACAATGTGCAGCTGGGTACATCCTTGAGCGAGAAGGCGGC
TGGGACAGTGTTCAAGACTGGATGGATGTACTCAGCGGAGGAGAGAAACAAGAATGGCGATGGCAAGAC
TGTTTTATCATAAACCCAGTTTGCCATTTGGATGAATGCACAAGTGGCGTCAGTGTGGATGTGGAGGA
CTACATTTACAGCCACTGTCGAAAGTTGGCATCACCCCTTTTACCGTTTACACAGAAAGTCCCTTTGG
AAGCACCAGGACTACTACCTGCATATGGATGGGAGAGGCAATTATGAATTCAAAAAGATCACAGAAGATA
CAGTTGAGTTCGGATCA

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR222851 representing NM_008991
Red=Cloning site Green=Tags(s)

MAAFSKYL TARNTSLAGAAFLLLCLLHKRRRALGLHGKKS GKPPLQNNEKEGKKERAVVDKVFLSRLSQI
 LKIMVPRTFCKETGYLLLIAMVLSRITYCDVWMIQNGTLIESGIIGRSSKDFKRYLFNFIAAMPLISLVN
 NFLKYGLNELKLCFRVRLTRYLYEEYLQAFYYKMGNLNRIANPDQLLTQDVEKFCNSVVDLYSNLSKP
 FLDIVLYIFKLTSAIGAQPASMMAYLLVSGFLFLTRLRRPIGKMTIMEQKYEYEGYRYVNSRLITNSEEIA
 FYNGNKREKQTIHSVFRKLVEHLHNFIFFRFSMGFIDSIIAKYVATVVGYLVVSRPFLDLAHPRLHSTH
 SELLEDYYQSGRMLLRMSQALGRIVLAGREMTRLAGFTARITELMQVLKDLNHGRYERTMVSQKEGIEG
 AQASPLVPGAGEIINTDNIKFDHVPLATPNGDILIQDLSFEVRSGANVLICGPNCGCKSSLFRVLGELW
 PLFGGRLTKPERGKLFYVQRPYMTLGLTRDQVIYPDGKEDQKKRGISDQVLKEYLDNVQLGHILEREGL
 WDSVQDWMVDVLSGGEKQRMAMARLFYHKPQFAILDECTSAVSVDVEDIYSHCRKVGITLFTVSHRSLW
 KHHEYLLHMDGRGNIEFKKITEDTVEFGS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9014_a10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_008991

ORF Size: 1977 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008991.3](#)

RefSeq Size: 3399 bp

RefSeq ORF: 1980 bp

Locus ID: 19299


UniProt ID: [P55096](#)

Cytogenetics: 3 52.94 cM

MW: 75.9 kDa

Gene Summary: The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the ALD subfamily, which is involved in peroxisomal import of fatty acids and/or fatty acyl-CoAs in the organelle. All known peroxisomal ABC transporters are half transporters which require a partner half transporter molecule to form a functional homodimeric or heterodimeric transporter. This peroxisomal membrane protein likely plays an important role in peroxisome biogenesis. Mutations have been associated with some forms of Zellweger syndrome, a heterogeneous group of peroxisome assembly disorders. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR222851