

## Product datasheet for **MR222851L4V**

### Abcd3 (NM\_008991) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Abcd3 (NM_008991) Mouse Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | Abcd3  |
| Synonyms:                 | AI313901; AU018866; AW146054; PMP; PMP68; PMP70; Pxm; Pxmp1  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_008991  |
| ORF Size:                 | 1977 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR222851).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_008991.2</a>  |
| RefSeq Size:              | 3399 bp  |
| RefSeq ORF:               | 1980 bp  |
| Locus ID:                 | 19299  |
| UniProt ID:               | <a href="#">P55096</a>   |
| Cytogenetics:             | 3 52.94 cM   |



[View online »](#)

**Gene Summary:**

The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the ALD subfamily, which is involved in peroxisomal import of fatty acids and/or fatty acyl-CoAs in the organelle. All known peroxisomal ABC transporters are half transporters which require a partner half transporter molecule to form a functional homodimeric or heterodimeric transporter. This peroxisomal membrane protein likely plays an important role in peroxisome biogenesis. Mutations have been associated with some forms of Zellweger syndrome, a heterogeneous group of peroxisome assembly disorders. [provided by RefSeq, Jul 2008]