

Product datasheet for **MR222852**

Rpgr (NM_001177954) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpgr (NM_001177954) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rpgr
Synonyms: Rd9; Rp3h
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR222852 representing NM_001177954
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCAAGAGGGTCGCGATGGGGTCCCAGGGGTAGGGCAGCACCTTAGGCTCAATCGAGTCGCCCTG
CTATCTTCCGAAGCAGGCACAGATTCCGTTTCGAGGCTTCGGCATGGCGGAATCTGAGTCACTGGTGCC
CGATACAGGTGCTGTGTTTACGTTTGGAAAACTAAATTTGCCGAAAATATTCCTAGCAAATCTGGTTT
AAAAATGACATACCCATATGTCTTTCATGTGGAGATGAACATACTGCTATTGTTACAGGAAATAATAAT
TGTACATGTTCCGCAGTAACAACCTGGGGTCAGTTAGGATTAGGATCAAAAGCTGCTATCATCAAGCCAAC
ATGTATCAAAGCTCTTAAGCCTGAGAAGGTGAACTTGCTGCCTGTGGAAGGAACCAACACCTTAGTTTCA
ACAGATACTGGTGGCGTATATGCAGCTGGTGGAAATAATGAAGGTCAACTGGGGCTTGGTGACACTGACG
ATAGAGACACCTTTCATCAAATGTCTTCTTTACACCTGCTGATACCATTAACAGCTCTCTGCTGGCGC
CAATACATCCGCTGCTCTTACTGAGGATGGAAAACCTTTTATGTGGGGTGACAATTCTGAAGGCGAGATT
GGTCTAGAAGATAAAAGTAATGTATGTATCCCTCATGAAGTGACTGTTGGAAAGCCAATTTCTGGATCT
CTTGTGGATATTACCATTCAGCTTTTGTAAAGAAGTAAGATGAACAGTTTTGATTCAATTTCCAAGTTATTA
CTGC

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR222852 representing NM_001177954
Red=Cloning site Green=Tags(s)

MPRGRSRWGSQGVGQHLRLNRVAPAIFFPKQAQIPFAGFGMAESES LVPDTGAVFTFGKTKFAENIPSKFWF
 KNDIPICLSCGDEHTAIVTGNKLYMFGSNNWGQLGLGSKAAIIKPTCIKALKPEKVKLAACGRNHTLVS
 TDTGGVYAAGGNNEQLGLGDTDDRDTFHQIVFFTPADTIKQLSAGANTS AALTEDGKLFMWGDNSEGQI
 GLEDKSNVCIPHEVTVGKPISWISCGYHSAFVTSKMNSFDSFSPSYCC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001177954
ORF Size: 774 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001177954.1, NP_001171425.1</u>
RefSeq Size:	2507 bp
RefSeq ORF:	777 bp
Locus ID:	19893
Cytogenetics:	X A1.1
MW:	28.3 kDa
Gene Summary:	<p>Could be a guanine-nucleotide releasing factor (By similarity). Plays a role in ciliogenesis (By similarity). Probably regulates cilia formation by regulating actin stress filaments and cell contractility (By similarity). May be involved in microtubule organization and regulation of transport in primary cilia (By similarity). Plays an important role in photoreceptor integrity. Isoform 5 may play a critical role in spermatogenesis and in intraflagellar transport processes. [UniProtKB/Swiss-Prot Function]</p>