

Product datasheet for **MR222856L4V**

Gabrb3 (NM_008071) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Gabrb3 (NM_008071) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Gabrb3 |
| Synonyms: | A230092K12Rik; AW049585; beta3; Cp1; Gabrb-3 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_008071 |
| ORF Size: | 1419 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR222856). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_008071.3 , NP_032097.1 |
| RefSeq Size: | 5568 bp |
| RefSeq ORF: | 1422 bp |
| Locus ID: | 14402 |
| UniProt ID: | P63080 |
| Cytogenetics: | 7 33.53 cM |



[View online »](#)

Gene Summary:

Ligand-gated chloride channel which is a component of the heteropentameric receptor for GABA, the major inhibitory neurotransmitter in the brain (PubMed:9108119). Plays an important role in the formation of functional inhibitory GABAergic synapses in addition to mediating synaptic inhibition as a GABA-gated ion channel (PubMed:27129275). The gamma2 subunit is necessary but not sufficient for a rapid formation of active synaptic contacts and the synaptogenic effect of this subunit is influenced by the type of alpha and beta subunits present in the receptor pentamer (PubMed:27129275). The alpha1/beta3/gamma2 receptor exhibits synaptogenic activity whereas the alpha2/beta3/gamma2 receptor shows very little or no synaptogenic activity (PubMed:27129275). Functions also as histamine receptor and mediates cellular responses to histamine (By similarity). Plays an important role in somatosensation and in the production of antinociception (PubMed:10670447).[UniProtKB/Swiss-Prot Function]