

## Product datasheet for **MR222857L3V**

### Chrnd (NM\_021600) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Chrnd (NM_021600) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Chrnd
Synonyms:	Ac; Achr; Achr-4; Acrd; L10076
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021600
ORF Size:	1563 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222857).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021600.2</a> , <a href="#">NP_067611.2</a>
RefSeq Size:	3003 bp
RefSeq ORF:	1563 bp
Locus ID:	11447
UniProt ID:	<a href="#">P02716</a>
Cytogenetics:	1 44.07 cM



[View online »](#)

**Gene Summary:**

This gene encodes the delta subunit of the muscle-derived nicotinic acetylcholine receptor, a pentameric neurotransmitter receptor and member of the ligand-gated ion channel superfamily. The delta subunit together with the alpha subunit forms the ligand-binding site. [provided by RefSeq, Nov 2012]