

## Product datasheet for **MR222894L4V**

### Lyar (NM\_025281) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Lyar (NM_025281) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lyar
Synonyms:	MLZ-264
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025281
ORF Size:	1167 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222894).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_025281.3</a> , <a href="#">NP_079557.2</a>
RefSeq Size:	1568 bp
RefSeq ORF:	1167 bp
Locus ID:	17089
UniProt ID:	<a href="#">Q08288</a>
Cytogenetics:	5 20.29 cM



[View online »](#)

**Gene Summary:**

Plays a role in the maintenance of the appropriate processing of 47S/45S pre-rRNA to 32S/30S pre-rRNAs and their subsequent processing to produce 18S and 28S rRNAs (By similarity). Also acts at the level of transcription regulation. Along with PRMT5, binds embryonic globin promoter (By similarity). Represses the expression of embryonic globin Hbb-y gene (PubMed:25092918). In neuroblastoma cells, may also repress the expression of oxidative stress genes, including CHAC1, HMOX1, SLC7A11, ULBP1 and that encoding the small nucleolar RNA SNORD41 (By similarity). Preferentially binds to a DNA motif containing 5'-GGTTAT-3' (By similarity). Stimulates phagocytosis of photoreceptor outer segments by retinal pigment epithelial cells (PubMed:25735755). Prevents NCL self-cleavage, maintaining a normal steady-state level of NCL protein in undifferentiated embryonic stem cells (ESCs), which in turn is essential for ESC self-renewal (PubMed:19489080).[UniProtKB/Swiss-Prot Function]