

Product datasheet for MR222934

Setdb1 (NM_001163642) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Setdb1 (NM_001163642) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Setdb1
Synonyms: AU022152; ESET; KMT1E; mKIAA0067
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR222934 representing NM_001163642
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGATCGCC

ATGTCCTCCCTCCCTGGGTGCATGAGTTGGCTGCAGCGCCAGCTGCAGCTGACTCTGCAGAGATTGCTG
AGCTGCAGCAGGCGGTGGTTGAAGAGCTGGGTATCTCTATGGAGGAACCTCGTCAGTACATTGATGAGGA
ACTGGAAGATGGACTGCATACAGCAGCGCAAGAAGCAGCTCGCAGAGCTGGAGACGTGGTACTACAG
AAAGAGTCTGAAGTGGCTTATGTTGATCGGCTGTTTATGATGATGCATCCAGGGAAGTGAAGTAACTGTGAGT
CTTTGGTGAAGGATTTCTACTCTAAGCTGGGACTACAGTATCATGACAGTAGCTCTGAGGATGAAGCTTC
CCGGCCACAGAGATCATTGAGATTCTGATGAAGATGATGATGTCTCAGTATTGATTCAGGTGATGCT
GGGAGCAGAACTCCAAAAGACCAGAAGCTTCGTGAAGCTATGGCTGCCTTAAGAAAATCAGCTCAAGATG
TCCAGAAGTTCATGGATGCTGTCAACAAGAAAAGCAGTTCTCAAGATCTACATAAAGGAACCTTGGGTCA
GGTGTCTGGGAACTGAGCAAAGATGGGGACCTGATAGTCAGCATGCGGATTCTGGGCAAGAAGAGGACT
AAGACATGGCACAAGGCACCCCTATTGCCATCCAGACTGTTGGGCTAGGAAAAAATACAAAGTGAAT
TTGACAACAAGGAAAGAGTCTGCTATCTGGGAACCATATTGCCATGATTACCACCCCTCCCGCTGACAA
GCTGTTTGTGGGAGTCGAGTGGTGGCCAAGTACAAGATGGAATCAGGTCTGGCTTATGCTGGCATT
GTAGCTGAGACCCCTAACGTCAAGAACAAGCTCAGATTTTTAATTTTTTTTATGATGATGGCTATGCTTCT
ATGTCAGTCAGTCAGAGCTTTATCCCATTTGCCGACCACTAAAAAAGACTTGGGAGGACATAGAAGATAG
CTCCTGCCGAGACTTCATAGAGGAATATCACTGCCTATCCAAACCGCCCAATGGTACTTCTCAAGAGT
GGGCAGCTTATCAAGACTGAGTGGGAAGGCACATGGTGGAAAGTCTCGAGTTGAAGAGGTGGATGGCAGCC
TAGTCAGGATCCTCTTTCTGGATGACAAAAGATGTGAGTGGATATATCGAGGCTCTACACGCTGGAACC
TATGTTTAGTATGAAGACATCCTCAGCCTCTGCAATGGAGAAGAAGCAAGGGGGCAACTCAGAACCCGT
CCTAATATGGGTGCTGTGAGGAGCAAAGTCTGTTGTTCAAGTATACACAGGATCTAACTGGTACTGGAA
TCCAGTTTAAGCCATGGAGCCCTACAGCCTATAGCTCCACCGCCCACTTCCATACCTCCTCTTTC
CCCCAAGCAGCTGACACTGACTTAGAAAAGCAACTTGACAATCACGGAACAAGTAGCCAAGAAGAGC
ACATCATCCGACCAGGATCTGTGGGCTCCGGCCATTCTCCCCTACTTCATCCACTCAGTGAATG



[View online >](#)

TGCTCTGCTGGGAAACTTGGGATAAACCAGACATATCGGTACACCTTTGGCCTCAGTAAACATCTACCCAGC
ATCTGCAGCCCCTCCAGTCCCTCCAGTCCCACCAGGGCCTCCAACCCCTCCAGGGCCTCCAGCTCCTCCA
GGGCTCTAGCTCCTCCAGCCTTCCATGGCATGTTAGAGCGGGCACCAGCTGAGCCCTCCTACCGAGCCC
CCATGGAGAAGCTTTTCTATTTACCTCATGTCTGCAGTTACACTTGTGTTGCCGGATCAGACCCATGAG
AAACGAACAGTATCGGGGCAAGAACCCTCTATTAGTCCACTTCTGTATGACTCCGGAGGATGACAGCA
CGGCGCAGAGTTAACCGCAAAATGGGCTTTCATGTAATCTATAAGACACCCTGTGGTCTCTGCCTTCGGA
CGATGCAGGAGATAGAGCGCTACCTTTTTGAGACTGGCTGTGACTTTCTGTTCTCCGGAGATGTTCTGTTT
GGATCCATATGTTCTTGTGACAGAAAGTTTCAACCCTTTAAGCCTTTTTACTATATTTTGGACATCACC
TATGGCAAGGAAGATGTTCCCTGTCTGTGTTAATGAGATTGACACAACCTCCCCACCCAGGTGGCCT
ACAGCAAGGAACGCATTCTGGCAAGGGTGTTCATTAAACACAGGCCCTGAATTTCTGGTTGGCTGTGA
CTGCAAGGATGGGTGTCGGGATAAATCCAAATGTGCCTGCCACCAGCTAACTATCCAGGCCACAGCCTGT
ACCCAGGGGGCCAAAGTCAACCCTAACTCTGGCTACCAGTATAAAAGACTAGAAGAGTGTCTGCCACAG
GGGTTTATGAGTGAACAAACGCTGCAATTGTGACCCAAACATGTGCACAATCGTTGGTGCAGCATGG
TCTGCAGGTTGACTACAGCTGTTAAGACACAGAACAAGGGCTGGGTATCCGCTGCTTGGATGATATT
GCCAAAGGCTCTTTGTCTGCATTTATGCAGGCAAAATCTGACAGATGACTTTGCAGACAAAGAAGGCC
TGAGATGGGTGATGAGTACTTTGCAAACTGGACCACATTGAAAGTGTGGAGAACTTCAAGGAAGGATA
TGAGAGTGTGTCCTCTGACAGCAGTGGGGTAGATATGAAGGACCAGGAAGATGGCAACAGC
GGTTCAGAGGACCTGAAGAATCCAATGATGACAGCTCTGATGATAACTTCTGTAAGGATGAGGACTTCA
GCACCAGTTCAGTGTGGCGTAGCTATGCTACCCGGAGGCAGACTCGGGGTCAAAGGAGAATGAATTGTC
TGAGATGACTTCCAAGGACTCCCGCCCCCAGACTCGGGCTCCACATGTTCTATCCCTTCTCAGTA
TCTGTAGGGGGTGCATCCACCTTCTCTGAAGACACCCAAGAACAAGGTGGCCTCGTGGTTGAGTT
GCAATAGTGTGAGTGAAGGTGGATTTGCTGACTCTGACAGCCGTTCTTCTTCAAGACTAGTGAAGGTGG
AGATGGCCGTGCTGGGGAGGCCGGGAGAGGCTGAAAGGGCCTTACCTCAGGATTGAGCTTCAAGGAT
GAAGGAGACAATAAGCAGCCTAAAAAGAGGACCCTGAGAACCAGAAACAAGATGCCAGTAGTTACTGAAG
GCTCTCAGAATCATGGACATAATCCTCCATGAAGTCTGAAGGGCTTCGCCGATCAGCTAGTAAATGTC
TGTGCTCCAGAGCCAGCGAGTTGTGACTTCTACTCAGTCAAACCCTGATGACATCCTGACTGTCCAGC
AGCACAGAGAGTGAAGGGGAAAGTGGAACCAGCCGAAAGCCACTGCTGGTCACTTACGCCACAGCTG
TTGATAGTGTGACATCCAGACCATCTTCTGGCTCTGACGGTGTGACTTTGAGGACAAGAAGAACTT
GTCAGGACCAACAAGCGCCAGGTGGCAGTAAATCAACCCGAGGCTTTGCTCTTAAATCAACCCATGGT
ATTGCCATTAATCAACCAACATGGCTTCCGTGGACAAGGGGAGAGTGCACCAGTTCGTAAGAACACAC
GCCAGTTCATGATGGTGAAGAGTCTTGCTACATCATTGATGCCAACTTGAAGGCAACCTAGGCCGCTA
CCTCAATCACAGTTGCAGCCCCAACCTGTTTGTCCAGAATGTGTTTGTGGATACCCATGATCTTCGCTTC
CCTTGGGTGGCCTTCTTTGCCAGCAAGAGAATCCGGGCTGGAACAGAACTCACTTGGGACTACAACACG
AAGTGGGCAGTGTGAAGGCAAGGAGCTGCTGTGCTGCTGTGGGGCCATTGAATGCAGAGGGAGACTTCT
T

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR222934 representing NM_001163642
 Red=Cloning site Green=Tags(s)


MSSLPGCMSLAAAPAAADS AEI AELQQAVVEELGISMEELRQYIDEELEKMDCIQQRKKQLAELETWV LQ
 KESEVAYVDR L FDDASREVTNCE SLVKDFYSKLG LQYHDSSEDEASRPT E I I E I P D E D D D V L S I D S G D A
 GSRTPK D Q K L R E A M A A L R K S A Q D V Q K F M D A V N K K S S S Q D L H K G T L G Q V S G E L S K D G D L I V S M R I L G K K R T
 K T W H K G T L I A I Q T V G L G K K Y K V K F D N K G K S L L S G N H I A Y D Y H P P A D K L F V G S R V V A K Y K D G N Q V W L Y A G I
 V A E T P N V K N K L R F L I F F D D G Y A S Y V T Q S E L Y P I C R P L K K T W E D I E D S S C R D F I E E Y I T A Y P N R P M V L L K S
 G Q L I K T E W E G T W W K S R V E E V D G S L V R I L F L D D K R C E W I Y R G S T R L E P M F S M K T S S A S A M E K K Q G G Q L R T R
 P N M G A V R S K G P V V Q T Q D L T G T G I Q F K P M E P L Q P I A P P A P L P I P P L S P Q A A D T D L E S Q L A Q S R K Q V A K K S
 T S F R P G S V G S G H S S P T S T L S E N V S A G L G I N Q T Y R S P L A S V T S T P A S A A P P V P P V P P G P P T P P G P P A P P
 G P L A P P A F H G M L E R A P A E P S Y R A P M E K L F Y L P H V C S Y T C L S R I R P M R N E Q Y R G K N P L L V P L L Y D F R R M T A
 R R R V N R K M G F H V I Y K T P C G L C L R T M Q E I E R Y L F E T G C D F L F L E M F C L D P Y V L V D R K F Q P F K P F Y I L D I T
 Y G K E D V P L S C V N E I D T T P P P Q V A Y S K E R I P G K G V F I N T G P E F L V G C D C K D G C R D K S K C A C H Q L T I Q A T A C
 T P G G Q V N P N S G Y Q Y K R L E E C L P T G V Y E C N K R C N C D P N M C T N R L V Q H G L Q V R L Q L F K T Q N K G W G I R C L D D I
 A K G S F V C I Y A G K I L T D D F A D K E G L E M G D E Y F A N L D H I E S V E N F K E G Y E S D V P T S S D S S G V M K D Q E D G N S
 G S E D P E E S N D S S D D N F C K D E D F S T S S V W R S Y A T R R Q T R G Q K E N E L S E M T S K D S R P P D L G P P H V P I P S S V
 S V G G C N P P S S E E T P K N K V A S W L S C N S V S E G G F A D S D S R S S F K T S E G G D G R A G G G R G E A E R A S T S G L S F K D
 E G D N K Q P K K E D P E N R N K M P V V T E G S Q N H G H N P P M K S E G L R R S A S K M S V L Q S Q R V V T S T Q S N P D D I L T L S S
 S T E S E G E S G T S R K P T A G H T S A T A V D S D D I Q T I S S G S D G D D F E D K K N L S G P T K R Q V A V K S T R G F A L K S T H G
 I A I K S T N M A S V D K G E S A P V R K N T R Q F Y D G E E S C Y I I D A K L E G N L G R Y L N H S C S P N L F V Q N V F V D T H D L R F
 P W V A F F A S K R I R A G T E L T W D Y N Y E V G S V E G K E L L C C C G A I E C R G R L L

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Cloning Scheme:



Plasmid Map:


ACCN: NM_001163642

ORF Size: 3921 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163642.1](#), [NP_001157114.1](#)

RefSeq Size: 4641 bp

RefSeq ORF: 3924 bp

Locus ID: 84505

Cytogenetics: 3 F2.1

MW: 145 kDa