

Product datasheet for **MR222950L4V**

Npnt (NM_033525) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Npnt (NM_033525) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Npnt
Synonyms:	1110009H02Rik; AA682063; AI314031; Nctn; POEM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_033525
ORF Size:	1734 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222950).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033525.2 , NP_277060.2
RefSeq Size:	3734 bp
RefSeq ORF:	1737 bp
Locus ID:	114249
UniProt ID:	Q91V88
Cytogenetics:	3 G3



[View online »](#)

Gene Summary:

Functional ligand of integrin alpha-8/beta-1 in kidney development. Regulates the expression of GDNF with integrin alpha-8/beta-1 which is essential for kidney development. May also play a role in the development and function of various tissues, regulating cell adhesion, spreading and survival through the binding of several integrins.[UniProtKB/Swiss-Prot Function]