

Product datasheet for **MR222955L3V**

Ebag9 (NM_019480) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ebag9 (NM_019480) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ebag9
Synonyms:	AI835379; Rcas1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019480
ORF Size:	642 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222955).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019480.4
RefSeq Size:	1606 bp
RefSeq ORF:	642 bp
Locus ID:	55960
UniProt ID:	Q9D0V7
Cytogenetics:	15 B3.2
Gene Summary:	May participate in suppression of cell proliferation and induces apoptotic cell death through activation of interleukin-1-beta converting enzyme (ICE)-like proteases.[UniProtKB/Swiss-Prot Function]



[View online »](#)